








## CONTROL PANEL

-  Power / Mute
-  Weight of Patient (kg)
-  Static Mode / Dynamic Mode
-  Low Pressure Alert LED
-  Pressure Monitoring LED (IPS)


## OPERATING INSTRUCTION

### Optimize Mattress Setting

- The default pressure setting is for Entrix 5 / Entrix MR mattress ( 30kg to 150kg )
- When using the Entrix Lux mattress, press  and  simultaneously for 3 seconds to switch between “default setting” and “Lux only setting”
- \* The system will reset to the default setting every time it is turned off by pressing the power key
- \* **Memory Recall Function** : Turning off the power unit by directly switching off the rocker switch will allow the power unit to memorize its setting for future use



### Power

- Plug the power cord into a wall socket
- Switch on the rocker switch on the side
- Press  to power on the unit, and then the green LED will light up
- \* The LED will change to yellow when the power cord is removed without switching off the rocker switch, or when there is an abnormal power outage

### Mute

- Press  to turn off the audible indicator



### Weight of Patient (kg)

- Press  to adjust the weight of patient
- **30kg to 150kg** : Entrix 5 / Entrix MR mattress
- **30kg to 180kg** : Entrix Lux mattress
- \* The weight scale is only indicative
- If the mattress seems too soft or too rigid, please adjust the comfort setting to conform to each patient’s requirements



## Therapy Mode Selection

Press  to switch between Static Mode and Dynamic Mode



## Static Mode

Static Mode is the therapy in which all air cells maintain constant pressure support and maximize patient's contact area to redistribute pressure



## Dynamic Mode

Dynamic Mode is the therapy in which air cells continuously alternate in an A-B-A-B pattern to relieve pressure and increase blood flow of the patient's tissue

- \* The 3 cells at the head are always excluded from the alternating
- \* Dynamic Mode will take 10 minutes to finish a complete cycle



## Low Pressure Alert LED

- An yellow LED will light up if the mattress has insufficient pressure
- \* The LED will stay on until the problem has been solved



## Pressure Monitoring LED (IPS)

- **IPS = Intelligent Pressure Sensing** The system monitors the mattress' pressure 24 hours
- The yellow LED indicates this function is in operation

| PROBLEM   | CONTROL PROCEDURE  | POSSIBLE SOLUTION  |
|---|--|--|
| The mattress is not inflating even though the power unit is working | <b>A</b> Verify that air flows freely across the tubes and the mattress manifold<br>Check for any cuts, blockages, or breaks | <b>A</b> If necessary, adjust the position of the tubes or manifold to prevent kinks or twists. Any cuts or rips in the air cells or air hoses, replace them |
|   | <b>B</b> Confirm that the quick coupling is properly connected to the air outlets of the power unit                          | <b>B</b> Make sure the quick coupling is firmly connected  |
|   | <b>C</b> Verify the CPR valve is correctly closed  | <b>C</b> Ensure the CPR valve is securely closed   |
| The patient sinks into the mattress                                 | <b>A</b> Check the weight setting on the power unit  | <b>A</b> Increase the patient's weight setting until the correct support pressure is achieved  |
|   | <b>B</b> Check for any abnormal air loss from the mattress   | <b>B</b> Replace the abnormally losing air components with an authentic replacement part   |
|   | <b>C</b> Check the air filter  | <b>C</b> Clean or replace the air filter   |
|   | <b>D</b> Verify that the CPR valve is correctly closed   | <b>D</b> Firmly close the CPR valve  |
| The power unit cannot power on                                      | <b>A</b> Verify that the power cord set plugs into the proper socket   | <b>A</b> Insert the power cord set of power unit into an appropriate socket and turn the power on  |
|   | <b>B</b> Verify that the power cord set is securely connected to the power unit  | <b>B</b> Insert the power cord set into the power unit and turn the power on   |
|   | <b>C</b> Verify that the power cord set is not damaged   | <b>C</b> Replace with a functioning power cord set   |
|   | <b>D</b> Verify that the fuses are not Blown   | <b>D</b> Contact the authorized distributor for technical service  |
|   | <b>E</b> The power unit is not responding to the control procedures listed above   | <b>E</b> Contact the authorized distributor for technical service  |