


POWERED BY 

What is CoolCore®

- Chemical Free material • UPF/Flammability rated • 4-Way Stretch nylon interlock with Elastane
- Patented combination of chemical-free blends that deliver 3 distinct functions:
Moisture Wicking to Moisture Transportation to Moisture Evaporation.
 - **CoolCore® Technology is exclusive to Stealth Products.**



Year Round Temperature Control

CoolCore fabrics wick moisture away from the body in every direction keeping you dry and comfortable no matter what season it is.

Breakthrough Technology

Research driven, patented fabrics from Coolcore are chemical-free and yield remarkable benefits. They cool better, last much longer and offer protection from the sun's harmful rays (UPF).

Stealth Products Cushion Specifications

Cushion Name	Weight Capacity	Weight of Cushion	Cushion Height	Width Range	Depth Range
Simplicity G	450 lbs	2.2 lbs	2"	16" - 24"	16" - 20"
Solution SPP	500 lbs	5.4 lbs	3"	10" - 24"	10" - 24"
Spectrum Foam SPP	650 lbs	3 lbs	3"	16" - 32"	16" - 28"
Spectrum Gel SPP	500 lbs	3.5 lbs	3"	14" - 24"	16" - 24"
Tru Comfort 2 SPP	450 lbs	4 lbs	4.5"	16" - 22"	14" - 24"
Essence	350 lbs	3.5 lbs	3.25"	16" - 20"	16" - 20"
Essence w/ Sm Wedge	350 lbs	4 lbs	3.5"	16" - 20"	16" - 20"
Essence w/ Med Wedge	350 lbs	4 lbs	3.5"	16" - 20"	16" - 20"
Essence Plus - Visco	350 lbs	4 lbs	4.5"	16" - 20"	16" - 20"
Essence Plus - HR	350 lbs	4 lbs	4.5"	16" - 20"	16" - 20"
Essence Plus - Gel	350 lbs	7 lbs	5"	16" - 20"	16" - 20"

*Please note: At the time of this printing all sizes may not be available. Please consult your local sales representative for current availabilities.

stealthproducts.com

Distributed by
 Pride Mobility Products Canada • 5096 South Service Road • Beamsville, Ontario L0R 1B3
 PH: 888-570-1113 • FAX: 866-514-1303 • www.pridemobility.com



Solutions For Life



POWERED BY 

stealthproducts.com

Stealth Products Cushions Color Guide

To make your cushion selection process a little easier, Stealth Products Cushions will all be marked with a color coded pvc patch. The color of the patch will identify the function and HCPC code of the cushion.



E2601/E2602



E2603/E2604



E2605/E2606



E2607/E2608



E2622/E2623



E2624/E2625



E2609

SIMPLICITY G



E2601/E2602

Features:

- Contoured, high-density molded foam base
- Coccyx and sacral relief
- Moisture resistant and easy to clean Stealth-Tek incontinence cover.
- Sizes available - 16" - 24"W x 16" - 20"D
- Mesh cover optional

SPECTRUM FOAM SPP



E2607/E2608

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Multi-layer, viscoelastic foam insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence, 4-Way Stretch Poly liner optional
- Sizes available - 16" - 32"W x 16" - 28"D
- Custom sizes available
- Optional Stealth-Tek, Mesh

SOLUTION SPP



E2607/E2608

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, aqueous gel insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- Sizes available - 10" - 24"W x 10" - 24"D
- Custom sizes available
- Optional Stealth-Tek or Mesh cover.

SPECTRUM GEL SPP



E2607/E2608

Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, lightweight fluid insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- 4-Way Stretch Poly liner optional
- Sizes available - 14" - 24"W x 16" - 24"D
- Custom sizes available
- Optional Stealth-Tek, Mesh

TRU-COMFORT 2 SPP

E2607/E2608



Convenient side pouch on both sides of the cushion with easy pull rubber loop zipper handles. Pouch is large enough to fit a larger format mobile phone.

Features:

- Deep contoured, high-density molded foam base
- CoolCore Technology cover standard.
- Easy pull rubber loop zipper handles
- Convenient side pouches on both sides of the cushion
- Sizes available - 16" - 22"W x 14" - 24"D
- Custom sizes available
- 4-Way Stretch Poly liner optional

ESSENCE SPP



E2607/E2608

Features:

- Available in Basic-Contour, Low-Contour and Medium-Contour configurations
- Lateral and pre-ischial shelf contoured, high-density specially formulated molded foam base
- Standard Visco topper for added comfort
- Coccyx and sacral relief
- CoolCore® technology cover standard.
- 4-Way Stretch Poly liner optional
- Sizes available - 16"W x 16"D - 20"W x 20"D
- Custom sizes available
- Optional Mesh outer cover



ESSENCE PLUS SPP



E2607/E2608

Features:

- Standard adjustable high profile wedge kit
- Lateral and pre-ischial shelf contoured, high-density specially formulated molded foam base
- Standard Visco topper for added comfort
- Coccyx and sacral relief
- CoolCore® technology cover standard.
- 4-Way Stretch Poly liner optional
- Sizes available - 16"W x 16"D - 20"W x 20"D
- Custom sizes available
- Adjustable contour wedge package provides better positioning customization.
- Topper options include, gel or HR foam
- Optional Mesh outer cover



Adjustable contour wedge package.