

Colibri

THE COLOURFUL ROUTE TO
EFFORTLESS, SAFE DRIVING



The Invacare Colibri is a stylish, simple and colourful micro lightweight scooter, designed for those who enjoy a totally independent lifestyle. Less than a minute is all that's needed to disassemble the Colibri for transportation.



Yes, you can.®

**Making Life's
Experiences Possible**

Colibri



Thanks to the unique patented Invacare LiteLock™ system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot. Plus, with its small footprint the Invacare Colibri offers greater manoeuvrability in tight, or limited spaces.

Additionally, you can choose from a range of trendy colour panels to match your style, swapping them according to your needs.

Choose where to charge your battery: directly on the Colibri or detached from it.

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.

Features



















Effortless dismantling can be done in less than a minute for transportation.



A soft padded seat provides additional comfort

Tech Specs

 Backrest width 465 mm	 Seat depth 400 mm	 Floor to seat height 370 - 420 mm	 Backrest height 385 mm	 Driving unit width 505 mm	 Total length 1010 mm	 Empty weight 49 kg (18 Ah)	 Max user weight 136 kg
 Battery capacity 2 x 18 Ah	 Turning diameter 2200 mm	 Max climable obstacle height 45 mm	 Ground clearance 50mm	 Max safe slope 6° / 10.5°	 Max driving distance Up to 16 km	 Max driving speed 8 km/h	 Motor capacity 200 W

Colour Options



* Warranty Terms and Conditions apply. Contact Invacare for full details.

FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz
Australia: Freephone 1800 460 460, email orders@invacare.com.au or visit www.invacare.com.au
Asia: Phone +61 2 8839 5330, email anzexportsales@invacare.com or visit www.invacare.com
Thailand: Phone +66 0 2821 5515 or visit www.invacare.com



Yes, you can.®