

# Introducing the **NEW GENIE LITE** T300

A super lightweight manual fold travel mobility scooter with **SUSPENSION**

The **New Genie Lite** is the latest advanced lightweight folding mobility scooter. It effortlessly blends style and comfort with great user performance and functionality. Comes with a full magnesium alloy frame for super lightweight, and rear suspension for added comfort. Easy lightweight lifting and loading into the Car!



Web: [www.portablemobility.com.au](http://www.portablemobility.com.au)

Email: [info@portablemobility.com.au](mailto:info@portablemobility.com.au)

PH: 1300 760 071

\*Exclusively imported and distributed by Portable Mobility

## Description

- Super lightweight magnesium frame
- Folds in 5 seconds!
- Big rear Suspension (for comfortable riding)
- 15km range
- 125 kg weight capacity
- LED Front lights
- Honeycomb tyres for a softer ride
- Only 16.85KG \*(without battery and armrests)
- Angle adjustable tiller
- Electromagnetic braking system
- Comfortable padded seat and backrest
- Lightweight flip up and removable armrests
- Armrests width adjustable
- Battery on board/off board charging
- Convenient high charging point position on side of displ
- Rolls along like a travel suitcase when folded

## Specifications

- Rear suspension
- Max speed 6km/h
- Range 15km
- Max weight capacity 125kg
- Folding size 376x450x690mm
- Unfold size 960 x 450 x 950 mm
- Max Gradient 12 degrees
- Turning radius 1.4m
- Driving system – Rear wheel drive via transaxle
- Honeycomb 7" Solid Rear Tyre
- Honeycomb 6" Solid Front Tyre
- Brushless motor
- Lithium Battery 24v 10AH
- Charger: 25.9V 2A (Output:25.9V 2A)
- Electromagnetic braking system
- Flip up and removable armrests
- Armrests width adjustable
- On board/off board battery charging
- 16.85KG weight \*(without battery and armrests)
- Colour options red and blue



\*this is a guide only, specifications might vary slightly

