Bandit

#1 PULL APART



Australia's most comfortable small scooter that easily disassembles for transport or storage. Standard features include a plush seat that swivels and slides so you can find a sweet ride position. Pneumatic tyres ensure a soft and quiet ride pushed along by a powerful 300W motor and backed by the quality guarantee of Top Gun Mobility's 5 Year Warranty.





Bandit Features



- Pull apart for easy transport
- Sliding seat adjustable to multiple height drivers
- 9 inch pneumatic tyres to cushion from rough driving surfaces
- All around suspension for driver comfort
- Best power to weight ratio in its class
- High visibility headlights for driver safety
- Swivel seat for easy mount and dismount

Specifications

Motor Power	300W
Max. Speed	6.8km/h
Operation Range	Up to 25km
Max. Climbing Grade	6-8° (depending on weight of user)
Turning Radius	110cm
Ground Clearance	7.5cm (+/- 2.5cm)
Max. Loading Weight	136kg
Overall Dimensions	108 x 50.5 x 87cm
Wheel Size - Front	9" pneumatic
Wheel Size - Rear	9" pneumatic
Weight	46.9kg
Battery	12V 20Ah x 2 (lead-acid)
Controller	PG 45A
Brake System	Electromagnetic
Free Wheel Mode	Yes
Suspension	Front & Rear
Lights	Headlights
Charger Size	24V 2A





NOTE: All specifications subject to change without notice.

- The information contained herein is correct at the time of publication;
 due to our commitment to constantly improve our products during
 development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery Condition and tyre pressure.

