

ottobock.

kidevo prime & kidevo mini.

Rigid frame wheelchairs for kids.

Kids are one of a kind – so are the products of the *kidevo* family. Each wheelchair is individually adapted to the child and can be further adapted in size when needed.

Thanks to the adjustable center of gravity, optimal driving behavior as well as an active and independent movement of the child can be encouraged from early on.

Due to the variety of options and the custom-made seat and back systems, each wheelchair becomes a unique product.

Highlights

- For an active and independent mobility of the child
- Easy to manoeuvre
- Custom-made seat and back systems
- Adjustable in three dimensions
- High variety of options
- Clever tilting function if required (t-version)



kidevo prime & kidevo mini. Rigid frame wheelchairs for kids.



Easy to use brakes

- All *kidevo* wheelchairs can be chosen with different child friendly brakes such as the drum brake. Combinations of brakes are also possible.



Custom-made seat and back system

- The seat and back system contains of a huge variety of different upholstered custom-made seat and back units.



Functional seat tilt

- The .t-versions have a functional and space-saving seat tilt without changing the center of gravity and thereby, maintaining an optimal driving behavior.



Functional back

- The back of the *kidevo prime* is optionally removable, flexible, adjustable in angle and further stabilized.



Robust foot plate

- For an optimal foot positioning, the *kidevo* wheelchairs come with a robust and flat one piece foot plate. In addition, there is a high variety of more positioning options to choose from.



Adjustable head supports

- The different types of head supports are adjustable in angle, height and depth.



High variety of rear wheels

- Each wheelchair model of the *kidevo* family provides three different rear wheel sizes in combination with various types of tires and handrims.



High variety of options

- Other options such as belts, therapy trays or grip ring knobs are available. In addition, individual spoke protectors can be chosen for further personalization.



Technical data

Feature	kidevo mini	kidevo prime
Max. user weight	50 kg	75 kg
Seat width range	16 – 32 cm	26 – 50 cm 28 – 50 cm*
Seat depth range	16 – 35 cm	20 – 50 cm
Tilt	-5 to + 40**	-10 to + 50**
Drive wheel sizes	20, 22 und 24"	22, 24 und 26"
Crash tested (acc. ISO 7176-19)	✓	✓
Weight	7,4 kg 9,0 kg*	17,4 kg 19,0 kg*

* for tilting variant (.t-version)

Frame colors glossy (RAL-colors)



Frame colors matt (RAL-colors)

