

Feather Light Strength,  
Built to endure!

# Airlite



**TopGun**  
**Mobility**  
*Mobility & Homecare Specialists*

Engineered with cutting edge combination of carbon fibre and magnesium. The dynamic duo ensures an incredibly lightweight frame without compromising on strength. The carbon fibre offers unmatched durability and resilience, while magnesium keeps the chair agile and easy to maneuver. Together, they create a powerchair that is as strong as it is light, providing both performance and reliability.

## Airlite Features

- Mesh fabric seat and backrest cushion, providing extra comfort.
- High strength carbon fibre frame.
- Lightweight and compact design for easy maneuverability.
- Suitable for indoors and on flat roads.
- Removable battery case for easy charging.



## Specifications

Weight (Without Battery)	12.5kg
Safe Working Load (SWL)	125kg
Operational Range	Up to 12km
Motor Power	24V 150W
Battery	24V/10Ah Lithium Battery (Optional 1PCS or 2PCS)
Charger	29.4V 2A MIC, AC PLUG
Turning Radius	90cm Minimum
Frame Material	Carbon Fibre + Maegnesium
Tyres Sizes Front	7" PU Solid Tyre
Tyres Sizes Rear	10" PU Solid Tyre
Seat (W x D x H)	44cm x 42 x 50cm
Seat Cushion Fabric	Mesh Fabric
Max Climbing Grade	8° (subject to user weight and surrounding conditions)
Dimensions (L x W x H)	104 x 58.5 x 90cm
Folded Dimensions (L x W x H)	86.5 x 30 x 76.5cm

- **Note:**
- All specifications subject to change without notice.
- The Information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.