

# PREMIER A3

**250KG SAFE LIFTING LOAD**  
**69CM SEAT WIDTH**



## POPULAR CUSTOM OPTIONS\*



Taller/shorter by adjusting seat height, seat depth, legrest, back height, arm height and giving an overall extended legrest 100cm+



Removable armrest



Waterfall back for gluteal shelf



Thicker/ softer footrest



Pressure Care Roho, equagel, flexigel



# PREMIER A3

250KG SAFE LIFTING LOAD

69CM SEAT WIDTH



LEG  
ELEVATION  
ABOVE HIP  
HEIGHT

A premium chair made of durable materials to hold a heavier and wider load for long periods, delivering long lasting comfort and independence to recline and stand unassisted by a carer. Three medical grade motors for maximum weight load and provides independent motors that allow the legs and back to move separately, including legs above hip height.

- ✓ Independent footrest & backrest
- ✓ Leg rest 85 mm higher than the body
- ✓ Heavy duty TRIPLE motor design
- ✓ Extra large size
- ✓ 2 cushions

## SPECIFICATIONS

Seat height : 50cm

Seat depth : 52cm

Seat width : 69cm

Arm height : 19cm

Back height : 72.5cm

## LUXURIOUS UPHOLSTERY

AMBASSADOR FABRICS RANGE  
– CHOICE OF 4 STANDARD COLOURS



CHARCOAL



MOCHA




PEBBLE GREY




STEEL BLUE

## CUSTOM FABRIC OPTIONS\* BARRIER



CARRFLEX  
STORM GREY



CARRFLEX  
BLACK



CARRFLEX  
NAVY

## VINYL



VOLCANO



METAL



ZODIAC



GLACIER



CHERRY



CEMENT



OASIS



MAHOGANY



CLAY



CASHMERE



VANILLA



YARRA

\*Custom options incur additional costs and lead time.