



VOGUE WALKER RANGE

aspire for..

Comfort

Stability

Style



VOGUE CARBON FIBRE SEAT WALKER



One of the lightest seat walkers on the market! Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super lightweight carbon fibre frame for a superior walking experience.



MODEL

WAF705350 | Vogue Carbon Fibre - Medium - Black
 WAF705360 | Vogue Carbon Fibre - Tall - Candy Red

FEATURES & BENEFITS shared by the Vogue Carbon Fibre and Lightweight Aluminium Seat Walkers



European Design

Ergonomic design offers personalised height adjustment, promoting correct posture



Compact

Turning circle for improved mobility in tight spaces



Simple

Folding mechanism locks for easy transportation and storage



Handbrakes

Lock for user safety and stability

	LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT** (kg)	BAG S.W.L (kg)
MED	680	600	695 - 920	240	460	550	225	8"	4.9	6.2	5
TALL	680	605	770 - 1040	240	460	620	225	8"	5.2	6.3	5

VOGUE LIGHTWEIGHT SEAT WALKER



This modern seat walker has a clean and crisp frame design, making it a standout when going for a walk or browsing the shops!



AVAILABLE COLOURS:



Red



Champagne

WEIGHS JUST OVER 5KG




Discreet Cable Management



Lightweight Aluminium Body
Improves manoeuvrability and portability

MODEL

- WAF705310RE | Vogue Lightweight Aluminium - Red
- WAF705310CH | Vogue Lightweight Aluminium - Champagne



Storage

Removable bag carries up to 5kg of essentials



Backrest

Deep seat and wide backrest for improved comfort and support



Bonus

Walking cane holder

1 YEAR LIMITED WARRANTY

AS/NZ Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT** (kg)	BAG S.W.L (kg)
655	620	785 - 920	240	465	540	220	8"	5.3	6.4	5

* Carbon Fibre Walker - Medium
** Includes accessories, bag and backstrap.

VOGUE ADVENTURE SEAT WALKER



Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super sized front wheels can manage any surface, perfect for the adventurous user!



**Luxury Look & Feel
Danish Designed**

AVAILABLE COLOURS:



Black

AutoLocate

Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate.



MODEL

WAF705450 | Vogue Adventure - Black

FEATURES & BENEFITS



European Design

Ergonomic design offer personalised height adjustment, promoting correct posture



Compact

Turning circle for improved mobility in tight spaces




Simple

Folding mechanism locks for easy transportation and storage



Handbrakes

Lock for user safety and stability



Storage

Removable bag carries up to 5kg of essentials



Bonus

Walking cane holder

1 YEAR LIMITED WARRANTY

AS/NZ Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)
770	695	900 - 1035	230	540	530	245	Front 11" Rear 9"	8.8	10	5

VOGUE FOREARM WALKER



For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



White / Black

Lightweight Aluminium Body
Improves portability and manoeuvrability




MODEL

WAF705400 | Vogue Forearm - White/Black

FEATURES & BENEFITS



European Design
Ergonomic design offer personalised height adjustment, promoting correct posture



Simple
Folding mechanism locks for easy transportation and storage



Bonus
Walking cane holder



Handbrakes
Lock for user safety and stability



Storage
Removable bag carries up to 5kg of essentials

1 YEAR LIMITED WARRANTY

AS/NZ Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	DEPTH (mm)	SEAT WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)
660	620	1020 - 1270	240	460	620	225	8"	9.1	10*	5

* Includes accessories, bag and backstrap.




Ideal for use within the home, the Indoor Walker is perfect for portable dining and in home mobility assistance.



AVAILABLE COLOURS:



Red / Black



Champagne / Black



Portable Dining
Safely collect and transport meals or hot drinks

On-cart Dining

Turn Around.
Slide Over Knee.
Lock & Enjoy.

MODEL

- WAF705100RE Vogue Indoor - Red/Black
- WAF705100CH Vogue Indoor - Champagne/Black

FEATURES & BENEFITS



European Design
Ergonomic design offer personalised height adjustment, promoting correct posture




Compact
Turning circle for improved mobility in tight spaces



Conveniently
Stands up when folded for flat storage



Handbrakes
Lock for user safety and stability



Storage
Removable bag carries up to 5kg of essentials



Bonus
Walking cane holder

1 YEAR LIMITED WARRANTY

AS/NZ Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

LENGTH (mm)	OVERALL WIDTH (mm)	HEIGHT (mm)	FOLDED WIDTH (mm)	WHEEL SIZE (inch)	MINIMUM FUNCTIONAL WEIGHT (kg)	TYPICAL FUNCTIONAL WEIGHT* (kg)	BAG S.W.L (kg)	TRAY S.W.L (kg)
650	540	860 - 970	270	6"	4.4	5.1	5	5

* Includes accessories, bag and tray.



Oxygen
Bottle Holder Bag
WAA695600



Night Light - LED Lamp
WAA695630



Single Hand Brake System* -
Compatible with Aspire Vogue Carbon Walkers
WAA695640

**Single Hand Brake System are not compatible with the Vogue Indoor Walker*



CRUTCHES, FRAMES & STICKS



SCOOTERS



OVERBED TABLES



BED STICKS & POLES



RAMPS



AIR MATTRESSES & FOAM MATTRESSES



COMMODOES, URINALS & BED PANS



OVER TOILET AIDS



SHOWER STOOLS & CHAIRS



RAISED TOILET SEATS



SHOWER COMMODOES



BATH BOARDS & BENCHES



ADJUSTABLE BEDS



POWER WHEELCHAIRS



PRESSURE CUSHIONS



LIFT & RECLINE CHAIRS



Distributed by:

ASBR21.43 130820

