



BANDIT-Z



Bandit-Z Features

- Accessories Included:
 - Phone Holder
 - Cup Holder
- Four colour panel options included - Purple, Orange, Green and Red.
- Zero turn - Turns on the spot for unbeatable manoeuvrability in tight spaces.
- Extra wide seat for superior comfort and support on every ride.



8.75kg



4.01kg



20.44kg



12.64kg



0.57kg

Specifications

Weight	45.8kg
Safe Working Load (SWL)	147Kkg
Motor Power	250W Brush Motor
Controller	Controller PG
Charger Size	2A
Battery	24V 12Ah (Lithium)
Tyre Size (Front)	8" Solid
Tyre Size (Rear)	8" Solid
Max. Speed	6km/h
Seat (W x D)	17.2 x 17cm
Lights	Front, Rear, & Side Lights
Turning	Zero Turning
Dimensions (L x W x H)	105.8 X 48.7 X 88.3cm



MOVE
MOBILITY