

# Frog Cane holder

Spring-loaded clamp

Large foam pad



A simple device designed to hold a walking stick upright and within reach while not in use.

The product attaches easily to any flat surface such as table, desk or shelf. Large foam pads protect the surface from scratching while ensuring a stable, safe grip.

The Frog has a spring-loaded opening which clamps to the walking stick once inserted. Suitable for most walking sticks with a maximum shaft diameter of 25 mm.

Colour may vary.

CODE	BL1706
SUITS STICK WIDTH	25 mm max
CONFORMS TO	ARTG176830
WARRANTY	1 year

