

## Sheepskin Rug Hospital Grade Natural

Cool in summer and warm in winter

Gold or Ivory

All natural sheepskin

Ideal for people required to spend long periods sitting or lying down, these hospital grade sheepskin rugs are prepared from specially selected sheepskins.

The wool characteristics maximise pressure relief for low risk pressure area care and maximum comfort.

With a wool length of approximately 25 mm, this luxurious rug has a very dense and springy wool fibre, which supports the body and allows air to flow freely helping to keep cool in summer and warm in winter. The rug can also absorb up to 30% of its weight in water without feeling wet.

Leather side is buffed to a lovely suede finish. Beautifully hand trimmed shape and packaged with full care instructions.

100% Australian Sheepskin. Machine washable at 60°C.

CODE	BL2023	BL2024
COLOUR	Ivory	Gold
WIDTH	680 mm	
LENGTH	1000 mm	
WEIGHT	1 kg	
CAPACITY	Unlimited	
WARRANTY	1 year	

