

## Silicone Fibre Heel Protect

Cotton/Polyester breathable cover provides maximum cool comfort



Silicone coating lubricates each fibre to reduce friction and shear force while preventing moisture penetration or retention



One size fits all, this product provides pressure relief and protection for heels and ankles.

The Heel Protect is constructed from silicone coated cylindrical fibres that retain their loft and promote airflow to the skin's surface for superior patient support and comfort. The silicone coating lubricates each fibre to reduce friction and shear force while preventing moisture penetration or retention, making it an ideal choice for geriatric or palliative care.

Covered with a breathable cotton/ polyester fabric, the heel protectors are fully machine washable with standard detergent.

- pressure relief for feet, heels and ankles whether in bed or wheelchair
- fluid resistant silicone fibres provides extra support and comfort
- Flame Retardant 16 CFR Part 1632
- Fully machine washable up to 70°

CODE	BL0637
HEIGHT	30 mm
WIDTH	100 mm
DEPTH	140 mm
CAPACITY	100 kg
COMFORMS TO	ARTG176835
WARRANTY	1 year

