



## DELUXE BEDSIDE COMMODE MAXI - CHAMPAGNE 550MM

Item No. KA500ZD55V03

Padded arms, seat and backrest for increased comfort . Hinged seat to allow access to toileting bowl . Front and side skirt to hide the bucket underneath . Available in 44cm, 50cm, 55cm 60cm and 65cm seat widths .

### DESCRIPTION

- Padded arms, seat and backrest for increased comfort .
- Hinged seat to allow access to toileting bowl .
- Front and side skirt to hide the bucket underneath .
- Available in 44cm, 50cm, 55cm 60cm and 65cm seat widths .

### SPECIFICATIONS

Seat Width	550 mm
Floor To Seat Height	500 mm
Product Weight	8 kg
Max User Weight	175 kg
Brand	K Care



## DELUXE BEDSIDE COMMUNE MAXI - CHAMPAGNE 600MM

Item No. KA500ZD60V03

Padded arms, seat and backrest for increased comfort . Hinged seat to allow access to toileting bowl . Front and side skirt to hide the bucket underneath . Available in 44cm, 50cm, 55cm 60cm and 65cm seat widths .

### DESCRIPTION

- Padded arms, seat and backrest for increased comfort .
- Hinged seat to allow access to toileting bowl .
- Front and side skirt to hide the bucket underneath .
- Available in 44cm, 50cm, 55cm 60cm and 65cm seat widths .

### SPECIFICATIONS

Seat Width	600 mm
Floor To Seat Height	500 mm
Product Weight	8 kg
Max User Weight	195 kg
Brand	K Care



## DELUXE BEDSIDE COMMUNE MAXI - CHAMPAGNE 650MM

Item No. KA500ZD65V03

Padded arms, seat and backrest for increased comfort . Hinged seat to allow access to toileting bowl . Front and side skirt to hide the bucket underneath . Available in 44cm, 50cm, 55cm 60cm and 65cm seat widths .

### DESCRIPTION

- Padded arms, seat and backrest for increased comfort .
- Hinged seat to allow access to toileting bowl .
- Front and side skirt to hide the bucket underneath .
- Available in 44cm, 50cm, 55cm 60cm and 65cm seat widths .

### SPECIFICATIONS

Seat Width	650 mm
Floor To Seat Height	500 mm
Product Weight	9.5 kg
Max User Weight	210 kg
Brand	K Care