



ETON

'Box in box' dual motor lift chair with headrest

Stonewash (Amber) 

Stonewash (Dune) 



Specifications

Weight Capacity	180 kg (400 lbs.)
Seat Width	53.3 cm (20.9")
Seat Depth	52 cm (20.4")
Seat-to-Floor	48.2 cm (18.9")
Chair Width	87.6 cm (34.4")
Chair Depth	95.3 cm (37.5")
Chair Height	106.7 cm (42")
Operation	USB Hand Wand

Description

Aside from its stylish appearance, the Eton also boasts the new integral 'box in box' mechanism. This fully enclosed feature means that nothing can inadvertently obstruct the chair while in its fully lifted position.

Quick Information

- Innovative 'box in box' technology
- Dual motor lift chair
- Independent headrest function
- Hand wand with USB connectivity
- Your choice of two earthy colours

THEOREM
POWER LIFT RECLINERS

MOVE
MOBILITY
movemobility.com.au

info@theoremconcepts.com • theoremconcepts.com