



<16 km
Range per Charge



136 kg
Weight Capacity



136.6 cm
Turning Radius



movemobility.com.au



Exceptional Portability

The Go Go® Carbon is the perfect blend of strength and portability. Built with a lightweight carbon fibre frame, this 4-wheel scooter offers exceptional durability while weighing just 16 kg without its (lithium) battery. It folds down in seconds for convenient transport and storage, and can be towed like luggage through airports and other environments. Powered by an NFC reader and digital display screen the Go Go® Carbon is designed for a lifestyle on the go.





Pride's EZ Turn Technology



Feather Touch Disassembly



Breathable Mesh Seat

RANGE PER CHARGE (MAX) ^{1,3}	Up to 16.7 km
MAXIMUM WEIGHT CAPACITY	136 kg
TURNING RADIUS ²	94 cm
MAXIMUM SPEED ¹	6.8 km/h
LENGTH & WIDTH ^{2,5}	L: 103.6 cm W: 49.4 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	Compact, foldable ⁸ W: 43 cm (17") D: 43 cm (17")
SEAT-TO-DECK HEIGHT ²	41.6 cm - 45.3 cm
SEAT-TO-GROUND HEIGHT ²	46.3 cm - 51.4 cm
GROUND CLEARANCE ^{2,7}	3.4 cm at transaxle

OVERALL WEIGHT (NOT INC BATTERY) ^{2,6}	37.6 kg
WEIGHT OF HEAVIEST PIECE ²	15.8 kg (front section)
SUSPENSION	N/A
3W/4W OPTION	4W
FRONT TYRES	2" x 7" solid
REAR TYRES	2.5" x 8" solid
BATTERY REQUIREMENTS ⁹	(2) 18AH
BATTERY CHARGER	Off-board, 2A
BATTERY WEIGHT ²	6 kg
COLOUR OPTIONS	Red Black Blue White

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
3. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity.
4. Does not account for anti-tip clearance.
5. Width of base (not including seat).
6. Includes standard seat weight.
7. Does not account for anti-tip clearance.
8. Solid durable plastic.
9. Deep cycle (AGM or gel cell type recommended).

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.

