



OVER TOILET FRAME WITH SEAT FLAP

Item No. KA410Z

Designed for users requiring additional seating support while toileting . White powder coated 22mm steel tube . Adjustable height and angle to suit individual users and uneven floors . Seat clips on and splash guard removes easily for quick assembly, cleaning or replacing .



DESCRIPTION

- Designed for users requiring additional seating support while toileting .
- White powder coated 22mm steel tube .
- Adjustable height and angle to suit individual users and uneven floors .
- Seat clips on and splash guard removes easily for quick assembly, cleaning or replacing .
- Anti-slip rubber tips on each leg for stability .
- Independently tested to AS/NZS 3973:2009 .
- Suitable for institutional and domestic use .

SPECIFICATIONS

Overall Width	510 mm
Overall Depth	600 mm
Min Overall Height	680 mm
Max Overall Height	805 mm
Seat Width	450 mm
Min Seat Height	480 mm
Max Seat Height	605 mm
Product Weight	5 kg
Max User Weight	110 kg
Brand	K Care