

Auto Folding Splitter Mobility Scooter

This lightweight remote-control auto folding travel mobility scooter incorporates advanced innovative design features. Mojo Lit quickly and easily splits into two parts! This allows light and easy lifting into the car boot!





1300 110 159

Web: https://www.thecarekiosk.com.au/

Email: thecarekiosk@gmail.com

*Imported by The Care Kiosk - dealer enquiries welcome



Description

- Strong lightweight frame
- Auto folding viaremote control and console
- Suspension for a comfortable ride
- · Splits in two in seconds for easy lightweight lifting into the car
- Telescopic and angle adjustabletiller
- Comfortable seat
- Flip up armrests, heightadjustable
- · Plane safe lithium batteries
- Speed control dial
- Front headlights
- Key switch
- Pull along when fully folded like a travel suitcase
- USB

Specifications

Folded Size: 51 cm x 48 cm x 68 cm
 Unfolded Size: 101 cm x 45 cm x 94 cm

Safe Gradient: 0-12
Turning Radius: 1.4 m
Max Capacity: 125kg
Max drive distance: 22 km
Max Speed: 0-6 km/hr
Motor: 24v 120w

Lithium Battery: Plane safe 2x 24v 7.5 ah = 24v 15ah
 Driving system: Rear wheel drive viatransaxle

Tyres: Honeycomb tyres
Suspension: Coil spring suspension
Charger: 29.4v 5A off board
Brake system: Electromagnetic brakes

Weight: 19 kg + 9 kg (frame + drive segment) w/o batteries

Colour option: White. red. blue



