



novis molift[®]
by Etac

Sit-to-Stand Solutions

Molift Raiser Pro

The ergonomic platform for a safe and active sit-to-stand transfer

Early mobilization and rehabilitation

Molift Raiser Pro allows the user to use their own muscle power during a transfer. The natural movement pattern stimulated and the ergonomic design makes the user feel safe during the transfer. Raiser Pro helps the user to rehabilitate in all health environments during every transfer so that mobilization can begin early, already in the intensive care unit.

The stability and maneuverability make Raiser Pro ideal in confined spaces, such as toilets.

Accessories



Safety strap – Raiser Strap

Provides support around the back or hips when standing up.

S/M Item no.....82508
L/XL Item no.....82509



Safety strap with sliding function – Raiser Strap+

When a user needs to be assisted, Raiser Strap+ can be used by one or two caregivers.

S/M Item no.....82528
L/XL Item no.....82529



Raiser Belt

Raiser Belt facilitates the work of caregivers who assist users up to a standing position. Can be used by one or two carers

S/M Item no.....82964 (82508 + 82963)
L/XL Item no.....82965 (82509 + 82963)

Features

- A versatile handle with a soft warm feeling, with an optimal grip for both the user and the caregiver.
- The handle has two different fixing heights for accessories.
- Soft padded legrest with side support, height adjustable with only one hand.
- A central brake mechanism with a one-foot touch that gives audible feedback when activated.
- Low step-on height of 35 mm.



| Product | Item no. | Material |
|-------------------|----------|--|
| Molift Raiser Pro | 16090204 | Powder coated aluminum, plastic, steel, chrome 3 |
| Dimensions | | |
| Height | 1200 mm | |
| Width | 540 mm | |
| Depth | 630 mm | |
| Weight | | |
| Total weight | 10,7 kg | |



For more information visit etac.com

Molift Quick Raiser 1 and Quick Raiser 2

Compact design and excellent manoeuvrability even in narrow spaces

Compact and easy

By virtue of its compact design and low weight, Molift Quick Raiser 1 and 2 are stable and easy to manoeuvre. The unique wheel position and near-perfect weight distribution of the base provides a small turning radius with easy manoeuvrability even in narrow spaces.

Quick Raiser 1 features an extra low base of just 7 cm, and with a narrow front of 34 cm, it fits well around toilets and under most beds, chairs and wheelchairs.

Quick Raiser 2 features a low base of just 10 cm, fits under most beds and chairs with the benefit of electrically adjustable legs to enable accessibility almost everywhere.

Natural pattern of movement

Both the Quick Raiser 1 and 2 have inclined columns that imitate the natural pattern of movement and fulfils the hoisting needs of users with some degree of weight-bearing capacity. Molift Quick Raiser 1 and 2 have a safe working load of 160 kg.



Molift Quick Raiser 1



Molift Quick Raiser 2

| | | | |
|------------------------------|-------------|--|-------------------------------|
| Molift Quick Raiser 1 | Item no. | Battery | NiMH 14,4 V - 2,6 Ah |
| U-shape arm | M04150 | Charge time | approx. 3 hours |
| Molift Quick Raiser 2 | | Hoisting speed 75 kg | 50 mm/sec |
| U-shape arm | M04400 | Rated performance | 40 lifts at 75 kg and 500 mmt |
| Chassis Weight QR1/QR2 | 30/36 kg | Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution. | |
| Dimensions QR1/QR2 | | | |
| Height | 1280 mm | | |
| Leg height | 70/100 mm | | |
| Chassis length | 955/1150 mm | | |
| Chassis width | 610 mm | | |
| External Width of wheel base | 1200 mm | | |
| Technical | | | |

Molift Quick Raiser 1 Features



Extra low base

The extra low base of 7 cm, and narrow front, 34 cm, fits well around toilets and under most beds, chairs and wheelchairs.



Multiple adjustable knee pads

The knee pads are height adjustable and can be turned around for better fitting and to give alternative support for the legs.



Battery indicator

Battery indicator on the hand control showing when charging is needed.



Unique wheel positioning

The lightweight construction combined with the unique wheel positioning provides exceptional manoeuvrability.

Molift Quick Raiser 2 Features



Electrically adjustable legs

A perfect alternative when the hoist is used together with chairs with no accessibility underneath, such as electric wheelchairs and recliners.



Adjustable knee pads

The knee pads are height adjustable and can be turned around for better fitting and to give alternative support for the legs.



Battery indicator


All Molift hoists have an integrated battery indicator that provides audio and visual notification.



Unique wheel positioning

The lightweight construction combined with the unique wheel positioning provides exceptional manoeuvrability.

Accessories



Molift RgoSling StandUp

Molift RgoSling StandUp is designed for a 2-point lifting arm. It is quickly and easily placed around the user. By virtue of the design it is ideal for use in toilet situations.

Features

- Perfect solution for hoisting and moving between different sitting positions.
- Large and firm support.
- Easy application.

| Molift RgoSling | Size | Waist (mm) | Item no. |
|------------------------------|------|------------|----------|
| StandUp 2-point | XS | 530-710 | 1720710 |
| | S | 640-850 | 1720720 |
| | M | 740-960 | 1720730 |
| | L | 930-1220 | 1720740 |
| | XL | 1150-1430 | 1720750 |
| StandUp 2-point with support | XS | 530-710 | 1720810 |
| | S | 640-850 | 1720820 |
| | M | 740-960 | 1720830 |
| | L | 930-1220 | 1720840 |
| | XL | 1150-1430 | 1720850 |

For more information regarding Active and StandUp lifting arms, see previous.. For complete variant matrix, contact your local distributor.

Material

Polyester



For more information visit etac.com

Molift Quick Raiser 205

An ergonomic sit to stand hoist with outstanding lifting capacity

Stable, safe and compact

Molift Quick Raiser 205 is a state of the art sit-to-stand hoist, with an electrically adjustable leg base. The hoist is comfortable, stable and safe with exceptional manoeuvrability. The hoisting of Molift Quick Raiser 205 imitates the natural pattern of movement and fulfils the transfer needs of users with some weight bearing capacity.

The unique wheel position, angled push bar and the extraordinary weight distribution of the base, provides a small turning radius. With the small castor option, it fits under most beds and chairs while the electrically adjustable legs enable accessibility almost everywhere. The fixed-leg option benefits a very small footprint.

Features

- Ergonomic design - the inclined column encourages a natural pattern of movement.
- Active 4-point Lifting arm - provides, together with RgoSling Active, extra support when hoisting from sit to stand.
- Easy to operate in moving and handling situations; toilet, sit-to-stand and stand-up training.
- Automatic Warning for periodic inspection and service.



Molift Quick Raiser 205

Model item no.

| Model | Battery | Adjustable leg |
|-----------------------------------|---------|----------------|
| Active 4-point, castor size Ø75mm | SLA | M29000 |

| | |
|------------------------|-------|
| Castor size Ø 75/50mm | 65 mm |
| Castor size Ø 100/75mm | 90 mm |

Accessories

Size

Item no.

| | | |
|-----------------------------------|-------------------------|---------|
| StandUp 2-point lifting arm | | 2910061 |
| StandUp slings | see matrix on next page | |
| Paediatric Footplate | | 2910065 |
| Set of inlays for footplate | | 2920151 |
| Leg support, bracket and cushions | Small | 2910059 |
| | Large | 2910056 |
| Set of cushions | Small | 2910076 |
| | Large | 2910077 |
| Lithium power pack and charger | | 2910041 |

Technical

| | |
|-----------------------|-------------------------|
| Battery SLA | SLA 24V 2,9 Ah |
| Battery Litium | LiFePO4 25,6V 2,4 Ah |
| Charge time | 6 hours |
| Rated performance | 100 hoists/75 kg/500 mm |
| Lifting Range | 600 mm |
| Turning radius | 1150 mm |
| Hoisting speed 75 kg | 35 mm/s |
| Hoisting speed 205 kg | 30 mm/s |

Weight

| | |
|-------------------|--------------|
| Total weight (FL) | 44 (42,5) kg |
|-------------------|--------------|

Dimensions

| | |
|-----------------------------|--------|
| External width of wheelbase | 620 mm |
| Internal width of wheelbase | 200 mm |
| Height of footboard: | |




Molift RgoSling Active

The ideal sling solution for active transfers with Molift Quick Raiser 205


Comfort and easy transfer

Molift Sit to Stand slings are padded and developed for active hoisting together with Molift Quick Raiser 205. The slings provides perfect lumbar support and are the ideal solution for hoisting and moving between different sitting positions. The slings have an efficient anti-slip lining to prevent the sling from sliding up. The waist belt provides additional support and stability.

Molift RgoSling Active is designed for a 4-point lifting arm and has sliding loops for the possibility to an optimal posture during the sit to stand movement. Molift RgoSling Active provides additional grip and support around the hips when standing. The wide back-support distributes weight and pressure evenly over a large surface, which makes it very comfortable for the user.



Molift RgoSling Active



Molift RgoSling StandUp

RgoSling Active and Quick Raiser 205 with 4-point lifting arm

Molift RgoSling StandUp is designed for a 2-point lifting arm. It is quickly and easily placed around the user. By virtue of the design it is ideal for use in toilet situations.

Features

- Perfect solution for hoisting and moving between different sitting positions.
- Self-adjusting sliding loops (RgoSling Active).
- Large and firm support.
- Easy application.

| Molift RgoSling | Size | Waist (mm) | Item no. |
|------------------------------|-----------|------------|----------|
| Active 4-point | XS | 500-720 | 1720910 |
| | S | 590-830 | 1720920 |
| | M | 730-950 | 1720930 |
| | L | 960-1210 | 1720940 |
| | XL | 1200-1390 | 1720950 |
| StandUp 2-point | XS | 530-710 | 1720710 |
| | S | 640-850 | 1720720 |
| | M | 740-960 | 1720730 |
| | L | 930-1220 | 1720740 |
| | XL | 1150-1430 | 1720750 |
| StandUp 2-point with support | XS | 530-710 | 1720810 |
| | S | 640-850 | 1720820 |
| | M | 740-960 | 1720830 |
| | L | 930-1220 | 1720840 |
| XL | 1150-1430 | 1720850 | |

Information

For more information regarding Active and StandUp lifting arms, see previous.. For complete variant matrix, contact your local distributor.

Material

Polyester



For more information visit etac.com



Proud distributor of Molift products in Australia.

novis ***molift***[®]

1300 738 885 / novis@novis.com.au
12/12 Mars Rd, Lane Cove West, NSW 2066

Scan QR code for more information
novis.com.au/molift



24022022F