

RIDE WITH POWER,
MOVE WITH STYLE!

Mustang Pro



MOVE
MOBILITY
movemobility.com.au

movemobility.com.au

Featuring a powerful 1400W motor for a smooth and efficient performance. The heated single seat ensures a warm, cosy ride no matter the weather. Stay in control with built in GPS and Geo Fencing technology, providing location tracking and safety boundaries. The one touch start button makes getting on the road effortless. while the SOS button adds an extra layer of security. A sleek digital display keeps all essential information at your fingertips, making every journey safe, smart, and comfortable.

Mustang Pro Features

- GEO Fencing
- GPS
- Heated Single Seat
- One Push Start Button
- Portable Belongings Hook
- Right, Left and Upper Atmosphere Light
- Rear View Mirror
- Lockable Rear Storage Box
- Free Wheeling Device
- Digital Display
- USB Charging
- Bluetooth (Music Capability)
- Radio
- Hazard Warning Lamp
- Mobile Phone Holder
- Drink Holder



Specifications

Weight	150kg
Safe Working Load (SWL)	220kg
Safe Working Load of Box	8kg
Colour Options	Ferrari Red, Lamborghini Yellow
Operational Range	Up to 70km
Motor	1400W
Controller	PG200A
Charger	PF10.OA
Charger Size	24V 10A
Battery	100AH X 2
Seat	Heated Single Seat
Tyre (Front)	14" Front Tyre
Tyre (Rear)	15" Rear Wide Tyre
Wheel Hub	Aluminium Alloy
Maximum Speed	10km/h
Maximum Climbing Grade	12° (subject to user weight and surrounding conditions)
Turning Radius	185cm
Freewheel Mode	Yes
Ground Clearance	140cm
Electric Brake	Warner 24V
LCD	6.2 inches
Dimensions (L x W x H)	165 x 77 x 128cm

- Note:
- All specifications subject to change without notice.
- The Information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.