



MOBILITY

movemobility.com.au

WHILL Model Fi



WHILL

The Perfect Travel Power Chair Has Arrived.

Lightweight, foldable, maneuverable. The Model Fi allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.



Foldable Design

Easily transport the Model Fi on buses, trains, and planes – or in the trunk of a car or taxi.*



When not in use, the Model Fi can be conveniently stored in a closet or out-of-the-way corner of the room.



Lightweight

Weighing only 52.9 lbs., the Model Fi is one of the lightest power chairs on the market today.**



Top Speed



Range



Obstacle Clearance



30.7" turning radius for tight turns

Connect Intelligently

The WHILL app allows you to check key device information and drive the Model Fi remotely.



Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model Fi will keep you safe while you ride.




Safely drives up or down 10° slopes.



Customization

Choose from five colors and a great selection of accessories to match your personal style and enhance your driving experience.

Color Options



- White ● Black ● Light Green
- Red ● Light Blue



Controller can be positioned on either the left or right side.

Accessories



Battery



Battery Charger




Lap Belt



Smart Key



Joystick



Basket



Cane Holder



Dust Cover



Side Bag

Specifications

Driving Range	Up to 12.4 Miles	Maximum Speed***	3.7 MPH
Foldable	Yes	Max Incline	10°
Turning Radius	30.7"	Obstacle Clearance	1.4"
Seat Sizes Available	18"	Width Between Armrests	Adjustable: 2 positions (16", 18")
Device Width	21.9" / 23.8"	Device Length	36.8" (18.3" when folded)
Device Weight	58.9 lbs	Charging Time	~ 5 hours
App Enabled	Yes, with remote control	Air Travel Approved	Yes
Max Weight Capacity	Up to 250 lbs	Swappable Arm Covers	5 Colors
Battery	Lithium-Ion (25.3V 10Ah)	Water Resistance	IPX4

***Acceleration and/or deceleration may change when driving on steep slopes.



WHILL Inc.

951 Mariners Island Blvd., Suite 300, San Mateo, CA 94404
 (844) 699-4455 us.info@whill.inc www.whill.inc



The WHILL Model Fi is not considered to be a medical device and has not been submitted to the Food and Drug Administration for review or clearance.

740-02222