

Etac Fresh bath board



Creating Possibilities

Fresh bath board has won a number of awards for its superb design and excellent functionality.

The original

Fresh has won many design awards and it has set the standard for bath boards in the industry. Its front edge, placed on the bathtub's outer side is wide to provide maximum support for transfers. The midsection is narrowed off to allow ease of access for intimate hygiene. Fresh inner side is designed to give place for soap and shampoo within comfortable reach.

Secure mounting

The edge 'stops' are easy to mount and secure by using the included plastic coin. Four friction linings make the bath board stay securely on the bathtub. Fresh is available in two lengths to ensure a perfect fit to various bath tubs.



Support handle
The orange handle design makes it easy to grip.



Fresh is the original that set the standard for bath boards in the industry with its curved design.

Description	Edge stop range	Weight	Item no.
Fresh Bath Board - 69 cm 27 1/4" with Orange Safety Handle	From 636 mm to 390 mm	2.2 kg 4.6 lbs	81600114
Fresh Bath Board - 74 cm 29 1/4" with Orange Safety Handle	From 686 mm to 432 mm	2.3 kg 5.1 lbs	81600024

Dimensions

Width: 36 – 27 – 30 cm 14 1/4 – 10 3/4 – 11 3/4"

Thickness: 4 cm 1 1/2"

Material

Overall: Polypropylene

Edge stops Polypropylene with TPE surface



reddot



max 85°C
max 185°F
for maximum of 3 minutes



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



Phone 1300 734 862
support@patienthandling.com.au