



Not available in fabrics shown.  
See Fabric Selections below.



## Pride Lift Chairs

*Experience the Luxury*

# LC-358XL Lift Chair

### Comfort and Style at a Great Value

Comfort and style complement one another perfectly in the LC-358XL from Pride® Lift Chairs. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, the LC-358XL Lift Chair delivers style and durability that will last for years to come. Plus, the LC-358XL offers numerous standard features including a 500 lbs. weight capacity, superb comfort and quiet and smooth performance, all at a great value.


#### SPECIFICATIONS

	LC-358XL
<b>Seat-To-Floor:</b>	51 cm (20")
<b>Seat Depth:</b>	56 cm (22")
<b>Seat Width:</b>	66 cm (26")
<b>Top of Back to Seat:</b>	74 cm (29")
<b>Weight Capacity:</b>	227 kg (500 lbs.)
<b>Back Style:</b>	Button-back

#### FABRIC SELECTION



**Burgundy**  
(vinyl)



**Latte**  
(fabric)

Note: Measurements may vary +/- 2.5 cm (1").  
Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.



**MOBILITY**

[movemobility.com.au](http://movemobility.com.au)

*Experience the Luxury*

**Pride**  
Mobility Products  
Australia Pty. Ltd.

**PRIDE**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST™ —

20-24 Apollo Drive  
Hallam, Victoria 3803  
[www.pridemobility.com.au](http://www.pridemobility.com.au)

Tel: 03 8770 9600  
Fax: 03 9703 2720

