

# RM410 Rotating Seat Shower Chair



The RM410 Rotating Seat Shower Chair is designed to improve safety, comfort, and independence during showering. Featuring a 360° rotating seat with secure lock positions, it allows easier transfers and safer movement in the bathroom.

## Key Features

- Oval shaped aluminium frame
- Drainage holes
- Removable armrests
- Non-slip stoppers
- 360 degree rotating seat
- Removable backrest
- Height adjustable extension legs

## Specifications

Weight Capacity (SWL)	136 kg
Product Weight	6 kg
Seat Width	51 cm
Seat Depth	33 cm
Seat Height	50 – 62 cm
Armrest Height	34 cm
Overall Height	85 – 97 cm
Overall Depth	55 cm
Overall Width	58 cm