

Safari



MOVE
MOBILITY



Top Gun Safari offers all the comfort of a luxury scooter with the added protection of our Clear Vision , Impact Resistant Weather Canopy. Powering the Safari is an unbeatable high torque motor and high capacity battery combination. When only the best will do! Safari – Ready to take on the elements.

TopGun
Mobility
WHEELS FOR LIFE

Safari Features

- High performance 1650W 4 pole, 2.2HP motor
- Clear vision impact resistant weather canopy
- Shock absorbing flexi tiller
- Adjustable tiller angle
- Up to 203kg weight capacity
- Deluxe captain seat with comfort lumbar support
- 14" front and rear sporty alloy rims
- Seat post suspension



Specifications

Motor Power	1650W
Max. Speed	10km/h
Operation Range	Up to 60km
Max. Climbing Grade	12° (depending on weight of user)
Turning Radius	220cm
Ground Clearance	21.5cm (+/- 2.5cm)
Max. Loading Weight	203kg
Overall Dimensions	153.5 x 81.6 x 172.5cm
Wheel Size - Front	14" pneumatic
Wheel Size - Rear	14" pneumatic
Weight	127kg
Battery	12V 120Ah x 2 (lead-acid)
Controller	PG S-200A
Brake System	Electromagnetic 10Nm
Free Wheel Mode	Yes
Suspension	Yes
Rear Mirrors	Left & Right
Charger Size	24V 8A



NOTE: All specifications subject to change without notice.

- The information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery Condition and tyre pressure.

TopGun
Mobility
WHEELS FOR LIFE