



SHOWER STOOL WITH ARMS AND PADDED SEAT

Item No. KA222ZAP

Designed for users preferring additional seating support while showering . White powder coated 22mm steel tube . Adjustable height and angle to suit individual users and uneven floors .

Padded seat clips on for easy assembly, change and cleaning . Anti-slip rubber tips on each leg



DESCRIPTION

- Designed for users preferring additional seating support while showering .
- White powder coated 22mm steel tube .
- Adjustable height and angle to suit individual users and uneven floors .
- Padded seat clips on for easy assembly, change and cleaning .
- Anti-slip rubber tips on each leg for stability .
- Independently tested to AS/NZS 3973:2009 .
- Suitable for institutional and domestic use .
- Available in 45cm and 50cm seat widths .

SPECIFICATIONS

Overall Width	510 mm
Overall Depth	380 mm

Seat Width	340 mm
Seat Depth	290 mm
Min Seat Height	510 mm
Max Seat Height	630 mm
Product Weight	2.5 kg
Max User Weight	110 kg
Brand	K Care
