

# ALPERTON

Dual motor lift lay-flat recliner with headrest and lumbar

Merino (Fresh) ●

Merino (Latte) ●



## Specifications

Weight Capacity	158.7 kg (350 lbs.)
Seat Width	43.1 cm (16.9")
Seat Depth	46.9 cm (18.4")
Seat-to-Floor	45.7 cm (17.9")
Chair Width	74.9 cm (29.4")
Chair Depth	93.9 cm (36.9")
Chair Height	99 cm (38.9")
Operation	USB Hand Wand

## Description

Smart and compact, the petite Alperton dual motor lift chair, with independent headrest and lumbar, combines class with function. Available in two stylish colours, the Alperton takes comfort to an exciting new level.

## Quick Information

- Independent headrest, legrest and lumbar function
- Footrest extension
- Steel seat box
- All close and lift button
- Hand wand with USB connectivity

**THEOREM**  
POWER LIFT RECLINERS

**MOVE**  
MOBILITY  
movemobility.com.au

info@theoremconcepts.com • theoremconcepts.com