

Bruno

MOVE
MOBILITY
movemobility.com.au

MOST VERSATILE LIGHT POWER CHAIR



Smart & smooth!
Safety smart
Bruno comes standard with intelligent braking to support its class leading 180kg load capacity. Coupling pneumatic tyres and soft ride suspension, Bruno gives you a smooth ride even over the rough stuff. Standard equipment includes an attendant control holder making this lightweight ride the perfect choice for a day out and about! Hop on and let Bruno take you for a ride.



TopGun
Mobility
WHEELS FOR LIFE

Bruno Features

- Powerful twin 250W motors
- World class quality
- 12" Pneumatic rear tyres for extra comfort
- Maximum load capacity 180kg
- 24V 10AH battery
- Drive range maximum 30km (with optional 2nd battery)
- Base weight only 24kg (without batteries)
- Intelligent braking system
- Optional accessories: Rear back pack, side bags, battery, headrest and leg extension holders

Specifications

Motor Power	250W x 2
Max. Speed	10km/h
Operation Range	Up to 15km
Max. Climbing Grade	9° (depending on weight of user)
Turning Radius	78cm
Ground Clearance	5.5cm (+/- 2.5cm)
Max. Load Capacity	180kg
Overall Dimensions	95.4 x 62.9 x 94.8cm
Folded Size	62.9 x 34.6 x 75.2cm
Wheel Size - Front	8" PU solid tyre
Wheel Size - Rear	12" pneumatic
Weight	24kg
Battery	12V 10Ah (Lithium)
Controller	40A Brushless Motor
Brake System	Intelligent Braking
Drive System	2 Motor, Rear Wheel
Suspension	Spring Absorber
Charger Size	24V 2A

