

World's Lightest Mobility Scooter!

Only 14.8kg!!!



MOBILITY

movemobility.com.au

GO *lite*



With its innovative cutting edge technology and attractive designs, this lightweight mobility scooter will give you years of dependable operation and will enhance the quality of your life by providing you with easy manoeuvrability and the mobility to experience an active daily lifestyle.

Go-Lite Features

- Unleash peak power with the 120W Brushless Motor.
- World's Lightest Mobility Scooter!
- Safe Working Load of 124kg.
- Easily maneuver through tight spaces with the nimble 8" tyres.
- Experience extended power with the 25.2V/20A Battery.
- Amazing operation range of up to 10km.
- Experience precision with the Electromagnetic Braking System.
- Incredible Folded Dimensions.

Specifications

Weight	16.3kg
Weight (Without Batteries)	14.8kg
Safe Working Load (SWL)	124kg
Colour Options	Cobalt Blue Red
Motor Power	120W Brushless Motor
Battery - Lithium	25.2V/10A or 25.2V/20A
Turning Radius	240cm
Maximum Speed	7km/h
Front & Rear Tyres	8" Solid Tyres
Ground Clearance	7.49cm
Maximum Climbing Grade (Subject to Weight & Conditions)	6-8°
Operational Range	Up To 10km (10Ah Battery)
Brake System	Electromagnetic Brake
Dimensions (L x W x H)	102 x 43 x 89cm
Folded Dimensions (L x W x H)	102 x 45.7 x 42cm

MOVE
MOBILITY
movemobility.com.au

- Note:
- All specifications subject to change without notice.
- The Information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

