

# Mustang

**MOVE**  
MOBILITY  
movemobility.com.au



Luxury, Power, Comfort. Along with an array of safety features the Mustang is the most comfortable ride to date. All round double sprung suspension, seat post suspension, luxury seat with adjustable lumbar support and multiple driver adjustment options. Completing the Mustang is the most powerful motor and battery combination on the Top Gun range. Go where you want in comfort and style. Mustang.

**TopGun**  
Mobility  
WHEELS FOR LIFE

## Mustang Features

- High performance 1650W 4 pole motor
- Safety turning sensor
- Shock absorbing flexi tiller
- Adjustable tiller angle
- USB charging port
- Up to 180kg weight capacity
- Luxury high back seat with lumbar support
- Rider comfort seat post suspension
- Dual tiller operation for ultimate ride comfort
- Safety park brake
- Digital dash display
- This is Australia's Rolls Royce in mobility scooters
- Front and Rear suspension, double sprung suspension

## Specifications

Motor Power	1650W
Max. Speed	10km/h
Operation Range	Up to 60km
Max. Climbing Grade	12° (depending on weight of user)
Turning Radius	200cm
Ground Clearance	17.5cm (+/- 2.5cm)
Max. Loading Weight	180kg
Overall Dimensions	146.1 x 69 x 131.7cm
Wheel Size - Front	13" pneumatic
Wheel Size - Rear	14" pneumatic
Weight	90kg
Battery	12V 100Ah x 2 (lead-acid)
Controller	PG S-200A
Brake System	Electromagnetic 10Nm
Free Wheel Mode	Yes
Suspension	Front & Rear
Rear Mirrors	Left & Right
Charger Size	24V 8A

