

TITAN



Our sleek new
TITAN scooter.
Built Tough For
The Adventurous!
This mobility
scooter will take
you to the next
level in outdoor
adventure.



MOVE
MOBILITY

Titan Features

- Experience unmatched power with the powerful 900W 4 pole Motor.
- Elevate your comfort with the 20" Captains Seat.
- Discover superior traction with 15" Rear tyres.
- Great for outdoor adventures.
- A massive operation range of up to 40km/h
- Ensure reliable stopping power with the Titan's manual handbrake.
- Free wheel mode for easy maneuverability.



Specifications

Weight	183kg
Safe Working Load	160kg
Colour Options	Black
Motor Power	900W (25:1) 4 Pole
Battery	12V 80Ah x 2pcs
Charger	8A - Off Board 120/240V, 50/60Hz
Controller	S-Drive 200 Amp
Turning Radius	135 - 200cm
Maximum Speed	10km/h
Front & Rear Tyres	350 x 120mm front, 15" rear
Seat	Sports 20" Captain Seat
Maximum Climbing Grade	10°
Driving Distance Range	25 - 40km
Brake System	Hand Brake - Manual (no locking)
Dimensions (L x W x H)	162 x 79 x 136cm



NOTE: All specifications subject to change without notice.

- The Information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

