

Windsor

POWERING YOUR INDEPENDENCE

MOVE
MOBILITY
movemobility.com.au



Top Gun
Mobility Windsor
toilet/commode
lift seat delivers a
stable, safe, easy to
use toilet transfer
solution that restores an
individual's ability to safely
toilet in the comfort of their home.

The 'sit to stand' motion is recognized as the most efficient way of getting to one's feet as it allows the application of the body's major muscle groups in the correct order. It mimics the rise motion of most "rise and recline" lift chairs so is a familiar movement for many people.

TopGun
Mobility
WHEELS FOR LIFE

movemobility.com.au

Windsor Features

- Maintain a safe body position whilst using the toilet
- Cleaner and easier functionality, with remote control included
- Comfortable, Stable lifting/lowering onto the toilet
- Reduced impact on joints
- Remote control for ease of use
- Removable portable bed pan for convenient bedside use
- Rear wheels for easy manoeuvrability (move from room to room if required)
- Increased personal confidence and prolongs your ability to live at home
- Battery powered device that enables users with restricted mobility to get on and off the toilet at the touch of a button
- SOS Alarm Button

Specifications

Dimensions (L x W x H)	60.5 x 52.5 x 79.5cm
Weight	22kg
Max. Loading Weight	200kg
Distance between handrails	50.6cm
Support times for battery full charge	150
Working Voltage	110-240V adaptive voltage
Shape	Support 33 degree lift and down



movemobility.com.au



NOTE: All specifications subject to change without notice.

- The information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.