

One of the
**World's
Lightest**

iLite



TopGun
Mobility
Mobility & Homecare Specialists

The iLite Rollator is one of Australia's Lightest Rollators! Built for strength without compromising on style or convenience. With an impressive safe working load of 150kg, it offers reliable support for everyday mobility. Lightweight yet durable, it features an easy to fold design, smooth rolling wheels, and ergonomic hand grips for comfortable use. Whether indoors or outdoors, the iLite Rollator is a trusted companion for confident movement.

iLite Features

- The detachable backrest adds comfort on demand and is easy to clean.
- The nylon seat, is durable lightweight and easy to clean - this seat provides extra daily comfort.
- Hard on the outside, soft and patterned inside - the wheels grip better for smoother, safer movement.
- The adjustable handles let you find your perfect height for better comfort and control.
- Keep your essential close - the built in storage bag is roomy, secure, and easy to access.

Specifications

Weight	5.59kg
Safe Working Load (SWL)	150kg
Colour Options	Blueberry Blue, Liquorice Black, Mango Orange, Nomad Grey
Height - Floor to Seat	55cm
Dimensions (W x D x H)	43.7 x 23 x 81-91cm
Warranty	1 Year Warranty



- Note:
- All specifications subject to change without notice.
- The Information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.