



Telescopic healthcare tables

The Telescopic Healthcare Tables feature a unique gas assisted telescopic height adjustment, giving them incredible versatility in a healthcare setting: switch from a coffee table to an overbed table, at the press of a button

FEATURES AT A GLANCE

- Pneumatic telescopic height adjustment activated with tabletop lever
- Low position suitable for use as coffee table or side table
- High position suitable for use as an overbed or chair table
- Extremely portable with four swivel castors and locking ability
- Choice of standard or round models
- Satin black powdercoated frame
- Easy hygiene and infection control with both models

Telescopic Healthcare Tables

COBALT

HEALTH



Features

- **Easy height adjustment:** pull the table top lever and apply slight downward or upward pressure to get the perfect height
- **Height adjustment range 440-830mm:** the only overbed table on the market which allows for such a height adjustment range
- **Versatility for small spaces:** the telescopic mechanism allows the table to be used as a side table or coffee table, or an overbed or over-chair table; this serves as two pieces of furniture in a resident's room, where space can be at a premium
- **Mobile, nimble and versatile:** high-quality swivel castors allow the table to be moved easily around the room or from room to room. All castors are locking
- **High quality materials:** this table features an extremely sturdy C-shape frame, with strong vertical column and a high quality MDF tabletop
- **Standard or round design:** the standard top has an extended base frame to allow it to fit around wide recliner chairs, the round top has a compact base frame to allow it to fit in to smaller spaces and fit around a day chair

Technical specifications

Standard Telescopic Table - model no.: A01-OBT-TELE-S

Table top: 790 x 390mm

Base frame: 890 x 435mm, internal dimensions: 810 x 400mm

Height range: 440 - 830mm

Shipping weight: 18kg

Carton dimensions: 910 x 520 x 300mm

Round Telescopic Table - model no.: A01-OBT-TELE-R

Table top: 700mm diameter (round)

Base frame: 646 x 435mm, internal dimensions: 570 x 400mm

Height range: 440 - 830mm

Shipping weight: 18kg

Carton dimensions: 910 x 700 x 300mm

Safe Working Load: 30kg (max load)

Cleaning: use a cloth slightly moistened with mild detergent - avoid using undiluted bleach or other surface cleaners

Warranty: Cobalt Health offers a 12 month warranty on our stand up workstations

cobalthealth.com.au

sales@cobalthealth.com.au
+61 2 8853 1100
Unit 34, 9 Salisbury Rd, Castle Hill
NSW 2154 Australia