

TITAN



Our sleek new
TITAN scooter.
Built Tough For
The Adventurous!

This mobility
scooter will take
you to the next
level in outdoor
adventure.

Built Tough
for the
Adventurous!

Zoom
Mobility
EASY LIVING

Titan Features

- Experience unmatched power with the powerful 900W 4 pole Motor.
- Elevate your comfort with the 20" Captains Seat.
- Discover superior traction with 15" Rear tyres.
- Great for outdoor adventures.
- A massive operation range of up to 40km/h
- Ensure reliable stopping power with the Titan's manual handbrake.
- Free wheel mode for easy maneuverability.



Specifications

Weight	183kg
Safe Working Load	160kg
Colour Options	Black
Motor Power	900W (25:1) 4 Pole
Battery	12V 80Ah x 2pcs
Charger	8A - Off Board 120/240V, 50/60Hz
Controller	S-Drive 200 Amp
Turning Radius	135 - 200cm
Maximum Speed	10km/h
Front & Rear Tyres	350 x 120mm front, 15" rear
Seat	Sports 20" Captain Seat
Maximum Climbing Grade	10°
Driving Distance Range	25 - 40km
Brake System	Hand Brake - Manual (no locking)
Dimensions (L x W x H)	162 x 79 x 136cm

NOTE: All specifications subject to change without notice.

- The Information contained herein is correct at the time of publication; due to our commitment to constantly improve our products during development, we reserve the right to alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

**Zoom
Mobility**
EASY LIVING

