



# There's an Airgo® rollator for your needs!



## Lightweight Rollator

An economical choice with all the basic features.

- Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Removable wire basket (4).
- 8" non-marring wheels are ideal for indoor/outdoor use (5).
- Reflective trim for added visibility (6).
- Weight capacity: 300 lb (136 kg).



Ergo-form handgrips

5' - 6' 2"	33" - 39"	22"	17"	24 1/2"	27 1/2"
Model: 700-940 Colour: Pacific Blue Net weight (without accessories): 16.2 lb					
AADL C061 Patent No.: 7,066,484					



## Adventure™ 6 and 8

Light and dependable!

- Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Soft basket doesn't need to be removed when folding (5).
- 6" or 8" non-marring wheels (6).
- Reflective trim for added visibility (7).
- Weight capacity: 300 lb (136 kg).



Ergo-form handgrips



Dual height frame technology

Standard 5' - 6' 2"	32" - 38"	21 1/2"	17"	23 1/2"	25 1/2"
Hemi 4'3" - 5' 5"	29" - 35"	18 1/2"	17"	23 1/2"	23 1/2"
Model: 700-941 (Adventure 6) Colour: Cranberry Net weight (without accessories): 14.1 lb					
Standard 5' - 6' 2"	33" - 39"	21 1/2"	17"	24 1/2"	27 1/2"
Hemi 4'3" - 5' 5"	30" - 36"	18 1/2"	17"	24 1/2"	25 1/2"
Model: 700-942 (Adventure 8) Colour: Titanium Net weight (without accessories): 16.2 lb					
ADP type 2 #MW2GM0040 (Adventure 6)		AADL C061			
ADP type 3 #MW3GM0060 (Adventure 8)		AADL C061			
Patent No.: 7,066,484					



## Ultra-Light 6

The lightest rollator available: a breeze to lift up and take anywhere!

- Easy-to activate brakes (2).
- Padded seat and backrest (3).
- Locks in open position for convenient storage (4).
- Soft basket doesn't need to be removed when folding (5).
- 6" non-marring wheels are ideal for indoor use (6).
- Reflective trim for added visibility (7).
- Weight capacity: 250 lb (113 kg).



Ergo-form handgrips



Multi-set height technology



Ultra-light frame

5' - 6' 2"	30" - 34"	19 1/4" - 23 1/4"	17"	25"	31"
Model: 700-950 Colour: Eggplant Net weight (without accessories): 10.5 lb					
ADP type 2 #MW2GM0050 AADL C061					



## NAVIGATOR

Full 2-in-1 transport chair and rollator.

- Easy-to activate brakes (2).
- Full contoured backrest simply flips for use as a transport chair or rollator (3).
- 8" non-marring wheels are ideal for indoor/outdoor use (4).
- Soft basket doesn't need to be removed when folding (5).
- Seatbelt that stores in the under seat storage bag.
- Weight capacity: 250 lb (113 kg).



Ergonomic armrests



Removable footrests



Folds ultra-compact

5' - 6' 2"	29" - 35"	21"	18.5"	26"	28"
Colours: Titanium Silver (793-710), Pacific Blue (793-720), Cranberry (793-740) Net weight (without accessories): 19.5 lb					
ADP type 3 MW3GM0015 AADL C061 (AADL C063 with 703-120 oxygen bag) Patent No.: 7,540,527					



## eXcursion™

The ultra-light side-fold with superior comfort & stability!

- Full contoured backrest provides superior back support (3).
- Fold by simply pulling the strap, the soft basket folds with the frame (5).
- A folding lock keeps the unit closed for secure lifting & storing.
- Signature ergo-form handgrips provide comfortable posture & absorb vibrations (1).
- Comfort-grip wheels positioned outside the frame for increased stability (6).
- NEW tool-free easy to adjust brake system.
- Weight capacity of 300 lb / XWD weight capacity of 400 lb.



Ergo-form handgrips



Ultra-light frame



Easy to fold



Collapsible basket

X23 (23" seat)	5.9' - 6'11"	33.5" - 39.5"	23"	17.5"	28.5"
Colours: Pacific Blue (700-917) Net weight (without accessories): 17.1 lb					
X20 (20" seat)	5' - 6' 2"	30.5" - 37"	20"	17.5"	28.5"
Colours: Pearl Black (700-920), Pacific Blue (700-921), Cranberry (700-922), Lime (700-923) Net weight (without accessories): 15.4 lb					
X18 (18" seat)	4' 3" - 5' 5"	29.5" - 35"	18"	17.5"	28.5"
Colours: Pearl Black (700-926), Cranberry (700-927) Net weight (without accessories): 15.4 lb					
XWD	5' - 6' 2"	30" - 36.5"	20"	19.5"	28.5"
Colours: Pacific Blue (700-928), Pearl Black (700-929) Net weight (without accessories): 16.4 lb					

ADP type 3 MW3GM0065 (eXcursion™ X18) AADL C061 ADP type 3 #MW3GM0065 (eXcursion™ X20) AADL C061  
ADP type 3 MW3GM0065 (eXcursion™ X23) AADL C061 ADP type 3 AA0000124 (eXcursion™ XWD) AADL C064



## COMFORT-PLUS™ and COMFORT-PLUS™ XWD

Simply the most comfortable rollator on the market.

- Easy-to activate brakes (2).
- Locks in open position for added safety and in folded position for convenient storage.
- Comfort-Grip indoor/outdoor 8" wheels (6).
- Reflective trim for added visibility (7).
- Comfort-Plus™ weight capacity: 300 lb (136 kg).
- Comfort-Plus XWD™ weight capacity: 400 lb (181 kg).



Ergo-form handgrips



Dual height frame technology



Ergo-seat & back



Extra-large shopping basket

Comfort-Plus™ Standard	5' - 6' 2"	33" - 39"	22.5"	17"	24.5"
Standard Hemi	4'3" - 5'5"	30" - 36"	19.5"	17"	24.5"
Colours: Iridescent Blue (700-946) Net weight (without accessories): 17.2 lb					
Comfort-Plus™ XWD	5' - 6' 2"	33" - 39"	22.5"	19"	26.5"
XWD Hemi	4'3" - 5'5"	30" - 36"	19.5"	19"	26.5"
Model: 700-949 Colour: Iridescent Blue Net weight (without accessories): 17.8 lb					
ADP type 3 #MW3GM0050 (Comfort-Plus™)		AADL C061		Patent No.: 7,066,484	
ADP type 3 #MW3GM0055 (Comfort-Plus™ XWD)		AADL C064		Patent No.: 7,066,484	



## duo™

The rollator that can transport you! Reversible backrest and flip-down footrest bar for use as a transport chair.

- Easy-to activate brakes (2).
- Locks in open position for added safety and in folded position for convenient storage (4).
- Soft basket doesn't need to be removed when folding (5).
- Comfort-Grip indoor/outdoor 8" wheels (6).
- Reflective trim for added visibility (7).
- Weight capacity: 250 lb (113 kg).
- Transport mode weight capacity: 195 lb (88 kg).



Ergo-form handgrips



Ergo-seat & back



Occasional use transport feature

Standard	5' - 6' 2"	32" - 38 1/2"	22"	17"	24 1/2"
Hemi	4'3" - 5' 5"	30" - 36 1/2"	20"	17"	24 1/2"
Model: 700-994 Colour: Midnight Blue Net weight (without accessories): 17.8 lb					
Model: 700-996 Colour: Midnight Blue Net weight (without accessories): 15.9 lb					

ADP type 3 #MW3GM0040 AADL C061 Patent No.: 7,066,484 7,628,411 ZL 200620147861.3 ZL 200920005465.0