



# MAXI BATHROOM RANGE

aspire for... Safety Comfort Support

SEAT WIDTH UP TO  
**775MM**




[www.aspirecare.com.au](http://www.aspirecare.com.au)  
1300 133 120




# aspire MAXI BATHROOM RANGE

## ASPIRE MAXI SHOWER STOOL



## ASPIRE MAXI OVER TOILET AID



## ASPIRE MAXI SHOWER CHAIR



### INTEGRATED FEATURES

1	Super heavy duty width adjustable frame to 775mm	4	Height adjustable legs with easy circlip adjustment points
2	Soft padding on seat and arms for postural support and comfort	5	Raised armrests support the user to sit down and stand up steadily
3	Reinforced aluminium frame accommodates bariatric users and is corrosion resistant	6	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability

Code	Description	Overall Height	Overall Width	Overall Depth	Seat Height	Sitting Width	Seat Depth
BTS118020	Aspire MAXI Shower Chair	670 - 1005mm	645 - 845mm	495mm	455 - 575mm	575 - 775mm	455mm
BTS119030	Aspire MAXI Shower Stool	670 - 790mm	645 - 845mm	495mm	455 - 575mm	575 - 775mm	455mm
BTT145020	Aspire MAXI Over Toilet Aid	670 - 790mm	645 - 845mm	495mm	455 - 575mm	575 - 775mm	455mm

www.aspirecare.com.au  
1300 133 120

