

OVERBED TABLES - SPRING/GAS/TILT

All feature our unique wide C-frame base, which provide stability and allow them to fit around lift chairs.



Shown with plastic formed top

Shown with HPL top

SPRING LIFT TABLE

- Spring assisted lift mechanism
- Large handscrew for ease of height adjustment
- Wide C-base to fit around most lift chairs
- Two lockable castors
- High-quality powdercoated steel frame
- Woodgrain table top in Beech or Cream
- Woodgrain matches bedside cabinet range
- Height adjustment range: 730 - 1000mm

Dimensions

Table top: 790 x 390 mm
Base frame external: 890 x 440 mm
Base frame internal: 810 x 400mm

Optional table tops

Plastic formed top, High Pressure Laminate

Weight

15kg

Order codes

Beech: A01-OBT-SPR-BCH
Cream: A01-OBT-SPR-CRM

GAS LIFT TABLE

- Gas assisted lift mechanism
- One-touch height adjustment
- Large handscrew to tighten and stabilise table top
- Wide C-base to fit around most lift chairs
- Two lockable castors
- High-quality powdercoated steel frame
- Woodgrain table top in Beech or Cream
- Woodgrain matches bedside cabinet range
- Height adjustment range: 730 - 1000mm

Dimensions

Table top: 790 x 390 mm
Base frame external: 890 x 440mm
Base frame internal: 810 x 400 mm

Optional table tops

Plastic formed top, High Pressure Laminate

Weight

15kg

Order codes

Beech: A01-OBT-GSL-BCH
Cream: A01-OBT-GSL-CRM

TILT TABLE

- Split tilt top with ledge for reading material or device
- Gas assisted lift mechanism
- Large handscrew to tighten and stabilise table top
- Wide C-base to fit around most lift chairs
- Two lockable castors
- High-quality powdercoated steel frame
- Woodgrain table top in Beech or Cream
- Woodgrain matches bedside cabinet range
- Height adjustment range: 730 - 1000mm

Dimensions

Table top: 790 x 390 mm
Base frame external: 890 x 440mm
Base frame internal: 810 x 400mm

Weight

17kg

Order codes

Beech: A01-TILT-GSL-UNI-BCH
Cream: A01-TILT-GSL-UNI-CRM