

ComfiMotion Collection

Luxurious adjustable bed collection for an enhanced sleeping experience



 aspire ComfiMotion®

ComfiMotion

collection

We believe that creating a supportive and luxurious sleep environment can significantly improve health and wellbeing.

Our ComfiMotion Collection is designed to blend seamlessly into the home environment, whilst improving comfort and independence.


Create Your ComfiMotion



Bed Function



Headboard Style



Mattress Support



Accessories

© Aidacare Pty Ltd 2024.

No part of this brochure may be copied, distributed, reproduced, scanned or stored on any electronic database, whether in whole or in part in any form or by any means.


Due care has been taken to ensure all the information contained in this brochure is correct at the time of printing. All measurements, pictures, colours and weight capacities are to be used as a guide only. We reserve the right to modify the design or appearance of any product displayed in this brochure without prior notice.

ComfiMotion is a registered trademark of Aidacare Pty Ltd.

Activ Care

Ultra Low Adjustable Bed

Function & Style Support



Supporting the modern, independent lifestyle.

The Activ Care bed combines the best of in-home care, functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

Features

- **Ultra Low Height:** Minimum platform height of 195mm reduces falls risk, improving user safety.
- **Customisable Safe Entry/Exit Height:** Programmable single button press ensures perfect bed height to maximise independence.
- **Underbed Clearance:** Provides clear access for patient lifters and overbed tables.
- **Battery Backup:** Provides confidence and maintains independence during power outages.
- **Night Assist:** Luminous underbed halo for safe entry and exit.
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.

Functions

Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks.

Head Tilt: Relaxes neck muscles, providing optimal support whilst reading or watching TV.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Auto-Regression in the Backrest: Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution.

Clinical Positioning: One touch raises back and legs to alleviate cardiac exertion, improving resting comfort.

Leg Raise: Promotes circulation, relieving tired, swollen, fatigued or painful limbs.

Memory Function: Allows the user to set and return to their preferred position at any time with one click of a button.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Seated Position: Assists with bed entry and exit and upright dining.

Trendelenburg (and Reverse): Can assist care givers when repositioning a user*.

Lockable Castors: Eight heavy duty dual castors.

Upholstery



Portobello Grey

Sizes (Length x Width x Height Range)

Long Single
2030 x 920 x 195 – 700mm **BEB046740**

King Single
2030 x 1070 x 195 – 700mm **BEB046750**

Long Double
2030 x 1370 x 195 – 700mm **BEB046760**

Queen
2030 x 1520 x 195 – 700mm **BEB046770**



▲ Companion Sleeping

See Configurations section for more companion options.



10 YEAR LIMITED WARRANTY

SWL
250kg
LS & KS

SWL
360kg
LD & Q


*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

Visit comfirmotion.com.au for Limited Warranty T&Cs

Care

Adjustable Bed

Function & Style Support



For comfort you deserve.

The Care adjustable bed blends seamlessly into the home without compromising on the therapeutic needs of the user or safety of partners or carers. Select from multiple position options designed to maximise independent mobility and improve sleep quality.

Features

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times.
- **Enclosed Frame Design:** Ensures optimal positioning when getting into bed and safely covers the moving base platform.
- **Safe-Lowering Technology:** Ensures that back, leg and Hi-Lo functions only lower under the force of gravity, reducing potential for accidents or injuries.
- **Box-on-Box Design:** Facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.



Functions

Hi-Lo Functionality: Assists with manual handling and helps to mitigate fall risks.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Clinical Positioning: Raises back and legs to alleviate cardiac exertion, improving resting comfort.

Lockable Castors: Four heavy duty castors.

Leg Raise: Promotes circulation, relieving tired, swollen, fatigued or painful limbs.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Seated Position: Assists with bed entry and exit and upright dining.


Trendelenburg (and Reverse): Can assist care givers when repositioning a user*.

Upholstery



*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

Variable Bed Height Options



▲ Companion Sleeping

Sizes (Length x Width x Height Range)

Long Single
2030 x 920 x 280 – 580mm¹ BEB046720

King Single
2030 x 1070 x 280 – 580mm¹ BEB046730

1. With Low Profile Bed Legs.



Visit comfirmotion.com.au for Limited Warranty T&Cs

Relax

Fixed Height Adjustable Bed

Function & Style Support



Supporting the modern, independent lifestyle.

The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.

Features

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- **Press & Hold Safety Handset:** Intuitive, large button interface enhances usability.
- **Box-on-Box Design:** Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- **German Engineered Motors:** For smooth, quiet adjustment.
- **Personalised Positioning:** Hand control allows users to find their perfect position to maximise comfort and relaxation.
- **USB Charge:** Connect and charge your devices conveniently.



Functions

Head Tilt: Relaxes neck muscles, providing optimal support whilst reading or watching TV.

Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.

Leg Raise: Promotes circulation, relieving tired, swollen and/or painful limbs.

Zero Gravity: Helps to relax the circulatory system and muscles of the lower back.

Memory Function: Allows the user to set and return to their preferred positions at any time with one click of a button.

Upholstery



Portobello Grey

Sizes (Length x Width x Height Range)

Long Single
2030 x 920 x 280 – 700mm **BEB046900**

King Single
2030 x 1070 x 280 – 700mm **BEB046910**

Long Double
2030 x 1370 x 280 – 700mm **BEB046920**

Long Queen
2030 x 1520 x 280 – 700mm **BEB046930**



▲ Companion Sleeping
See Configurations section for more companion options.



10 YEAR LIMITED WARRANTY



Visit comfirmotion.com.au for Limited Warranty T&Cs


Mattresses

Engineered for Motion

Our ComfiMotion Pocket Spring and Memory Foam mattresses are designed with the flexibility to follow the contouring profile of an adjustable bed.

Pocket Spring

The ComfiMotion Pocket Spring mattress offers zoned support springs that respond and contour to the individual body shape to improve comfort, maximising support and sleep quality.




Transitional Layered Pocket Springs arranged in weight activated support zones.



Integrated Side-Walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.



Luxuriously Soft Pillowtop Layer provides a plush immersive feeling complimented by the responsive pocket springs.



▲ Pocket Spring Mattress

▲ Classic Memory Foam Mattress

Classic & Pressure Care

The ComfiMotion Classic and Pressure Care are memory foam mattresses that contour to the user's body, providing personalised comfort and transitional pressure relief. Integrated side-walls increase roll resistance at the bed edge and assist users, partners or carers with extra rigidity to ease bed entry and exit.



Contoured Immersion: Specially layered memory foam provides contoured immersion for support and pressure relief.



Superior Breathability: Up to 8x more breathable than conventional memory foams.



Anti-microbial: Environmentally safe treatment is added to the foam which reduces odours and kills 99.9% of bacteria.



Pressure Care Incontinence Cover
Breathable, 4-way stretch waterproof cover to reduce shear, friction and maintain skin micro-climate.



Integrated Side-Walls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit.

Pocket Spring

(Length x Width x Height)

Pocket Spring Mattress – Long Single
2030 x 920 x 240mm [BEM469800](#)

Pocket Spring Mattress – King Single
2030 x 1070 x 240mm [BEM469810](#)

Pocket Spring Mattress – Long Double
2030 x 1370 x 240mm [BEM469820](#)

Pocket Spring Mattress – Queen
2030 x 1520 x 240mm [BEM469830](#)



Visit comfimotion.com.au for Limited Warranty T&Cs

Memory Foam - Classic

(Length x Width x Height)

Long Single Classic
2030 x 920 x 200mm [BEM469910](#)

King Single Classic
2030 x 1070 x 200mm [BEM469920](#)

Long Double Classic
2030 x 1370 x 250mm [BEM469922](#)

Queen Classic
2030 x 1520 x 250mm [BEM469924](#)

Memory Foam - Pressure Care



(Length x Width x Height)

Long Single Pressure Care
2030 x 920 x 200mm [BEM469930](#)

King Single Pressure Care
2030 x 1070 x 200mm [BEM469940](#)

Long Double Pressure Care
2030 x 1370 x 250mm [BEM469942](#)

Queen Pressure Care
2030 x 1520 x 250mm [BEM469944](#)



Visit comfimotion.com.au for Limited Warranty T&Cs

Clinical Considerations

Accessibility and Inclusivity: Unexposed bed surround with padded upholstery removes gaps between beds in co-sleeping arrangements and promotes safe entry and exit.

User Independence: Height adjustability improves user mobility and assists carers with manual handling.


Flexible Positioning: Caters to various user requirements such as respiratory conditions, snoring, sleep apnoea, pressure redistribution, joint or muscle pain, lymphoedema and swelling.

Vision and Hand Function: Handset with a clear and intuitive button display, assists users with reduced vision or hand function.

Falls Prevention: Integrated bed lighting and a minimum platform height of 195mm (Activ Care) and 280mm (Care and Relax) helps to reduce the likelihood and severity of falls in the home.

User Safety: Battery backup provides peace of mind in the event of a power outage (Standard on Activ Care, optional accessory on Care).


Headboards



Allure

Stylish freestanding headboard
Subtle square padded design

Long Single **BEA830120**
King Single **BEA830130**
Long Double **BEA830132**
Queen **BEA830134**
King **BEA830136**



Chesterfield

Stylish freestanding headboard
Chesterfield padded design

Long Single **BEA830100**
King Single **BEA830110**
Long Double **BEA830112**
Queen **BEA830114**
King **BEA830116**

Pillows



Breeze Cooling Gel Memory Foam

- Adaptive and supportive sleep surface
 - Helps regulate sleep surface temperature
 - Accommodates side and back sleepers
- BEP050100**



Plush Memory Foam and Fibrefill

- Familiar sleep surface
 - Improved head and neck support
 - Accommodates side and back sleepers
- BEP050110**


Accessories

Function & Style Support




Crook Handle Bed Stick

Care Bed **BEA830200**
Activ Care Bed **BEA830211**
Provides multiple grip points to assist repositioning in bed, and safe entry and exit.




Extendable Side Rail

Activ Care Bed **BEA830205**
Provides a telescoping bed rail support for users at risk of falls.



Bed Rail

Care Bed **BEA830210**
Activ Care Bed **BEA830212**
Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility.



Split King Upgrade Kit

Activ Care Bed LS & KS only **BEA830230**
Ensures simultaneous Hi-Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.



Simplified Handset

Activ Care Bed **BEA830214**
For clients with reduced cognitive ability.



Partner Handset

Activ Care Bed **BEA830213**
Allows double and queen beds to be operated from either side of the bed.



Standard




Large

Large Accessibility Castors

Care Bed only (4 x required) **MWT451080**

Larger castors raise the bed frame to 155mm, facilitating overbed table and patient lifter access.



Battery Backup

Care Bed
*Standard on Activ Care Bed **BEA830220**

Provides peace of mind in the event of a power outage.


Collection Summary

Beds



Care Adjustable Bed – (LS) _____	BEB046720	Activ Care Adjustable Bed – (LS) _____	BEB046740
Care Adjustable Bed – (KS) _____	BEB046730	Activ Care Adjustable Bed – (KS) _____	BEB046750
Relax Adjustable Bed – (LS) _____	BEB046900	Activ Care Adjustable Bed – (LD) _____	BEB046760
Relax Adjustable Bed – (KS) _____	BEB046910	Activ Care Adjustable Bed – (Q) _____	BEB046770
Relax Adjustable Bed – (LD) _____	BEB046920	Fixed Companion Bed – (LS) _____	BEB046800
Relax Adjustable Bed – (Q) _____	BEB046930		

Mattresses




Pocket Spring Mattress – (LS) _____	BEM469800	Memory Foam Mattress – Pressure Care (LS) _____	BEM469930
Pocket Spring Mattress – (KS) _____	BEM469810	Memory Foam Mattress – Pressure Care (KS) _____	BEM469940
Pocket Spring Mattress – (LD) _____	BEM469820	Memory Foam Mattress – Pressure Care (LD) _____	BEM469942
Pocket Spring Mattress – (Q) _____	BEM469830	Memory Foam Mattress – Pressure Care (Q) _____	BEM469944
Memory Foam Mattress – Classic (LS) _____	BEM469910		
Memory Foam Mattress – Classic (KS) _____	BEM469920		
Memory Foam Mattress – Classic (LD) _____	BEM469922		
Memory Foam Mattress – Classic (Q) _____	BEM469924		

Headboards



Allure Headboard – (LS) _____	BEA830120	Chesterfield Headboard – (LS) _____	BEA830100
Allure Headboard – (KS) _____	BEA830130	Chesterfield Headboard – (KS) _____	BEA830110
Allure Headboard – (LD) _____	BEA830132	Chesterfield Headboard – (LD) _____	BEA830112
Allure Headboard – (Q) _____	BEA830134	Chesterfield Headboard – (Q) _____	BEA830114
Allure Headboard – (K) _____	BEA830136	Chesterfield Headboard – (K) _____	BEA830116

Pillows



Breeze Pillow _____	BEP050100	Plush Pillow _____	BEP050110
---------------------	---------------------------	--------------------	---------------------------

Accessories




Crook Handle Bed Stick – Care Bed _____	BEA830200	Large Accessibility Castors – Care Bed _____	MWT451080
Crook Handle Bed Stick – Activ Care Bed _____	BEA830211	Battery Backup – Care Bed _____	BEA830220
Bed Rail – Care Bed _____	BEA830210	Split King Upgrade Kit – Activ Care Bed _____	BEA830230
Bed Rail – Activ Care Bed _____	BEA830212	Simplified Handset – Activ Care Bed _____	BEA830214
Extendable Bed Rail _____	BEA830205	Partner Handset – Activ Care Bed _____	BEA830213

Configurations for Companions

Whichever ComfiMotion package you select, rest assured there is a companion solution available to suit those wishing to sleep alongside any ComfiMotion Adjustable bed.

What to Order


KING
2030 x 1840



KING PLUS
2030 x 1990



SUPERKING
2030 x 2140



Combine any two singles below:

Care Long Single _____	BEB046720
Active Care Long Single _____	BEB046740
Fixed Long Single _____	BEB046800
Relax Long Single _____	BEB046900


Combine any single with king single below:

Relax Long Single _____	BEB046900
Care Long Single _____	BEB046720
Active Care Long Single _____	BEB046740
Fixed Long Single _____	BEB046800
Relax King Single _____	BEB046910
Care King Single _____	BEB046730
Active Care King Single _____	BEB046750

Combine any two king singles below:


Care King Single _____	BEB046730
Active Care King Single _____	BEB046750
Relax King Single _____	BEB046910

Long Double and Queen Beds can not be paired with other beds in the range.




Dual Adjustable Beds

Designed to sit flush when positioned together to enable close monitoring and enhanced care. This configuration offers the flexibility of independent adjustability, without compromising on the sense of inclusivity and intimacy between partners.



Split King Upgrade Kit

Ensures simultaneous Hi-Lo platform movement whilst still allowing individual positioning requirements when beds are in partner configuration.



Fixed Companion Bed

Matching single bed without any adjustability designed to sit flush when positioned together. It facilitates a side-by-side sleeping arrangement at a fixed height of 370mm, without compromising on the users care and assistance requirements.

Split King Upgrade Kit [BEA830230](#)

Compatible with active care bed

Fixed Companion Bed – (LS) 2030 x 920 x 370mm [BEB046800](#)

LS – Long Single

KS – King Single

LD – Long Double

Q – Queen

K – King

AIDACARE PTY LTD

Building 3A, 1 Moorebank Avenue,
Moorebank NSW 2170
Australia

1300 133 120 | aidacare.com.au

AIDACARE LIMITED

Arcadia House, Maritime Walk,
Ocean Village, Southampton SO14 3TL
United Kingdom

01622 541235 | aidacare.co.uk

MADE IN PRC

ASBR28.95 160124

