


### Lift Seat Powered Sit-To-Stand and Vertical Toilet Lifts

- Lowers you from a standing position
- · Raises you back up from a sitting position
- Safely use the toilet on your own
- An easy-to-use, two-button controller
- Operate over your own toilet
- Reduces the risk of falls & injuries









## Our range of powered toilet lifts cater for all range of disabilities and have the following key features:

Works over any toilet and guaranteed\* to fit your toilet space (\*lift seat home swl 23st)

Standard and bariatric models all fully height adjustable

No installation, ground fixing or assembly required

Rechargeable battery and mains power options

Simple two button hand control - elongated toilet opening

Can be used bedside with our drop-in commode bucket

Works over Geberit<sup>1</sup> and Clos-o-Mat<sup>2</sup> wash dry toilets (all LIFT SEATs work over Geberit and the LIFT SEAT HOME PLUS and VERTICA work over Clos-o-Mat)

Raises and lowers smoothly - Safe simple and easy to use



## The Lift Seat Range of Powered Toilet Lifts



#### Lift Seat Home

**Lift Seat Home** (SWL 23ST/145KGS)

**Lift Seat Home Plus** (SWL 32ST/205KGS)

NATURAL SIT-TO-STAND MOVEMENT



#### Lift Seat Vertica

**Lift Seat Vertica** (SWL 23ST/145KGS)

**Lift Seat Vertica Plus** (SWL 32ST/205KGS)

VERTICAL UP AND DOWN MOVEMENT

#### **Dimensions**

#### Lift Seat Home


# 510mm / 20 in —

460mm / 18 in-395mm / 15.5 in

-535mm / 21 in **-**




#### Lift Seat Home Plus



#### Lift Seat Vertica

.....





The maximum height of all our products are fully adjustable with our levelling feet



0800 043 9358



liftseat.co.uk



info@liftseat.co.uk



@liftseat